CELL-DYN 26 Plus Control

CONTROL L N H

INTENDED USE

CELL-DYN 26 Plus Tri-Level Control is an assayed hematology control for evaluating the accuracy and precision of hematology instruments that provide a white blood cell differential.

SUMMARY AND PRINCIPLE

It is an established laboratory procedure to use stable quality control material to monitor the performance of diagnostic tests. CELL-DYN 26 Plus Control is composed of stable materials that provide a means of verifying the accuracy and precision of CELL-DYN hematology results. It is handled and run in the same manner as patient specimens, and it is available in three levels representing low, normal, and high results.

REAGENT

CELL-DYN 26 Plus Control is an in vitro diagnostic product that may contain any or all of the following: stabilized human or mammalian red blood cells, human, mammalian or simulated white blood cells and a platelet component in a preservative medium.

WARNINGS AND PRECAUTIONS

IVD

For In Vitro Diagnostic Use

CAUTION: This product contains human sourced and/or potentially infectious components. Refer to the **REAGENT** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens¹. Biosafety Level 2² or other appropriate biosafety practices^{3,4} should be used for materials that contain or are suspected of containing infectious agents.

- All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For additional information, refer to the appropriate CELL-DYN Operator's Manual.

INSTRUCTIONS FOR USE

Refer to the CELL-DYN 26 Plus Control Assay Sheet for mixing and handling instructions.

PROCEDURE

Refer to the quality control procedures provided in the appropriate CELL-DYN Operator's Manual.

STORAGE AND STABILITY

CELL-DYN 26 Plus Control should be tightly capped and stored at 2–10°C (36–50°F). Protect containers from overheating and freezing. Unopened containers are stable until the expiration date indicated. Once opened, containers can be used only for the number of days stated on the assay sheet, provided that they are handled properly: avoid unnecessary cycles of warming and cooling, prolonged exposure to ambient room temperature, or vigorous mixing, all of which may damage the control. In addition, the volume in the tube must meet or exceed the minimum sample volume stated in the Operator's Manual.

Trending in the MCV or RDW parameters over the product shelf life is inherent to hematology control products. This characteristic does not indicate product instability. Ranges and limits for these parameters may need to be adjusted.

INDICATIONS OF DETERIORATION

CELL-DYN 26 Plus Control, after mixing, should be similar in appearance to fresh whole blood. In unmixing, refrigerated containers, the supernatant fluid may appear cloudy and reddish. Other discoloration of the supernatant fluid, or marked hemolysis, may indicate deterioration. Inability to recover expected values may also indicate product deterioration. Rough handling, freezing, overheating, and contamination are frequent causes of product damage. Incomplete mixing or instrument malfunction may also cause unacceptable results. **Do not use the product if deterioration is suspected;** contact Abbott Customer Service.

Customer Service: Contact your local representative or find country specific contact information on www.abbottdiagnostics.com.

PRODUCT LIMITATIONS

Proper storage and use of this product as described in this document is required for optimal performance. This product is not suitable for a manual differential analysis of white blood cells. Values are not transferable to systems not listed on the assay sheet.

CELL-DYN Enzymatic Cleaner Concentrate solution must be thoroughly removed from the analyzer Wash Block and/or Probe before running control products.

EXPECTED VALUES AND THEIR DERIVATION

Refer to the table of values on the enclosed assay sheet for expected results. The assay values are obtained from replicate testing on CELL-DYN Systems. The mean range is an estimate of observed interlaboratory variation due to reagent differences, maintenance, calibration, and operating technique. The Systems are calibrated as follows:

- using whole blood according to CLSI and ICSH-recommended manual, reference methods^{5,6,7,8,9,10}
- using a commercially available CELL-DYN calibrator
- using the appropriate, reliably calibrated CELL-DYN hematology analyzer

PERFORMANCE CHARACTERISTICS

The assigned values are presented as a mean and a recovery range. The mean assay values are derived from repetitive testing on several instruments operated and maintained according to the manufacturer's instructions; they do not necessarily apply to a single instrument. The recovery ranges are intended to reflect inter-laboratory and inter-instrument variability; thus, they are wider than the ±2 SD QC range for one instrument. Always perform quality control according to good laboratory practice, laboratory director's requirements, and any regulatory or accreditation requirements.

CONTROL PERFORMANCE INQUIRY

If you experience control recovery problems, please have the following information available when you call Abbott Customer Service:

- Expiration dates, lot numbers and supporting data for control(s) in question.
- Previous control lot numbers and data for these previous lots.
- Data from current precision study (N=10) using fresh whole blood specimens and performed according to instrument operator's manual.
- Data from Quality control reports as well as data from your last instrument calibration.
- Instrument maintenance history.

BIBLIOGRAPHY

- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.
- The Expert Panel in Cytometry on the ICSH, J.M. England, et al.: *Recommended Methods for the Visual Determination of White Cell and Platelet Counts*. World Health Organization, WHO/ LAB/88.3, 1988.
- The Expert Panel in Cytometry of the ICSH, J.M. England, et al., *Reference Method for the Enumeration of Erythrocytes and Leucocytes*. Clin. Lab. Haemat. 16:131-138, 1994.
- Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard - Third Edition*; CLSI Document H15-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2000.
- Procedure for Determining Packed Cell Volume by the Microhematocrit Method; Approved Standard - Third Edition*; CLSI Document H07-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2000.
- The Expert Panel in Cytometry on the ICSH, J.M. England, et al.: *The Assignment of Values to Fresh Blood Used for Calibrating Automated Blood Cell Counters*. Clin. Lab. Haemat. 10:203-212, 1988.
- Validation, Verification, and Quality Assurance of Automated Hematology Analyzers*; Approved Standard - Second Edition; CLSI Document H26-A2, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2010.

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1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.
5. The Expert Panel in Cytometry on the ICSH, J.M. England, et al.: *Recommended Methods for the Visual Determination of White Cell and Platelet Counts*. World Health Organization, WHO/ LAB/88.3, 1988.
6. The Expert Panel in Cytometry on the ICSH, J.M. England, et al.: *Reference Method for the Enumeration of Erythrocytes and Leucocytes*. Clin. Lab. Haemat. 16:131-138, 1994.
7. *Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard - Third Edition*; CLSI Document H15-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2000.
8. *Procedure for Determining Packed Cell Volume by the Microhematocrit Method; Approved Standard - Third Edition*; CLSI Document H07-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2000.
9. The Expert Panel in Cytometry on the ICSH, J.M. England, et al.: *The Assignment of Values to Fresh Blood Used for Calibrating Automated Blood Cell Counters*. Clin. Lab. Haemat. 10:203-212, 1988.
10. *Validation, Verification, and Quality Assurance of Automated Hematology Analyzers*; Approved Standard - Second Edition; CLSI Document H26-A2, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2010.

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BIBLIOGRAFIA

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA; 2005.
5. The Expert Panel in Cytometry on the ICSH, J.M. England, et al.: *Recommended Methods for the Visual Determination of White Cell and Platelet Counts*. World Health Organization, WHO/ LAB/88.3, 1988.
6. The Expert Panel in Cytometry on the ICSH, J.M. England, et al.: *Reference Method for the Enumeration of Erythrocytes and Leucocytes*. Clin. Lab. Haemat. 16:131-138, 1994.
7. *Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard - Third Edition*; CLSI Document H15-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2000.
8. *Procedure for Determining Packed Cell Volume by the Microhematocrit Method; Approved Standard - Third Edition*; CLSI Document H07-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2000.
9. The Expert Panel in Cytometry on the ICSH, J.M. England, et al.: *The Assignment of Values to Fresh Blood Used for Calibrating Automated Blood Cell Counters*. Clin. Lab. Haemat. 10:203-212, 1988.
10. *Validation, Verification, and Quality Assurance of Automated Hematology Analyzers*; Approved Standard - Second Edition; CLSI Document H26-A2, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2010.

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