CELL-DYN 29 Plus Control (with Retic)

CONTROL L N H

INTENDED USE

CELL-DYN 29 Plus Control (with Retic) is a Whole Blood hematology quality control material used to monitor results obtained on CELL-DYN hematology systems.

SUMMARY AND PRINCIPLE

It is an established laboratory procedure to use stable quality control material to monitor the performance of diagnostic tests. CELL-DYN 29 Plus Control (with Retic) is composed of stable materials that provide a means to verify the accuracy and precision of results for hemogram, WBC differential, and reticulocyte parameters; it is handled and run in the same manner as patient specimens and is available in three levels representing low hemogram/high reticulocyte, normal hemogram/intermediate reticulocyte, and high hemogram/low reticulocyte results.

REAGENT

CELL-DYN 29 Plus Control (with Retic) is an *in vitro* diagnostic product that may contain any or all of the following: stabilized human or mammalian red blood cells, human, mammalian or simulated white blood cells, and a platelet component in a preservative medium.

WARNINGS AND PRECAUTIONS

IVD

For *In Vitro* Diagnostic Use

CAUTION: This product contains human sourced and/or potentially infectious components. Refer to the **REAGENT** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens'. Biosafety Level 2^o or other appropriate biosafety practices^{5,6,7,8,9,10} should be used for materials that contain or are suspected of containing infectious agents.

- Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. For detailed information visit [www.corelaboratory.abbott](#)→Customer Portal→Technical Library→Other Reference Documents→Hematology Aids.

Safety Data Sheets are available at [www.abbottdiagnostics.com](#) or contact your local representative.

For additional information, refer to the CELL-DYN Operator's Manual.

INSTRUCTIONS FOR USE

Refer to CELL-DYN 29 Plus Control (with Retic) Assay Sheet for mixing instructions.

PROCEDURE

Refer to the quality control procedures provided in the appropriate CELL-DYN Operator's Manual.

STORAGE AND STABILITY

Protect vials from overheating and freezing. CELL-DYN 29 Plus Control (with Retic) is stable through the expiration date when stored at 2° to 10°C (36-50°F). After opening, **CELL-DYN 29 Plus Control (with Retic)** is stable throughout the open-vial dating, as indicated on the Assay Sheet, when stored at 2° to 10°C. Once opened, vials can be used only for the number of days stated on the assay sheet, provided that they are handled properly: avoid unnecessary cycles of warming and cooling, prolonged exposure to ambient temperature, or vigorous mixing, all of which may damage the control. In addition, the volume in the vial must meet or exceed the minimum sample volume stated in the Operator's Manual.

Trending in the MCV or RDW parameters over the product shelf life is inherent to hematology control products. This characteristic does not indicate product instability. Ranges and limits for these parameters may need to be adjusted.

INDICATIONS OF DETERIORATION

CELL-DYN 29 Plus Control (with Retic), after mixing, should be similar in appearance to fresh whole blood. In unmixed, refrigerated containers, the supernatant fluid may appear cloudy and reddish. Other discoloration of the supernatant fluid, or marked hemolysis, may indicate deterioration. Inability to obtain expected values may indicate product deterioration. Rough handling, freezing, overheating, and contamination are frequent causes of product damage. Incomplete mixing or instrument malfunction may also cause unacceptable results. **Do not use the product if deterioration is suspected;** contact Abbott Customer Service.

Customer Service: Contact your local representative or find country specific contact information on [www.abbottdiagnostics.com](#)

PRODUCT LIMITATIONS

Proper storage and use of this product as described in this document is required for optimal performance. The components of this product are not suitable for microscopic differential analysis or spun microhematocrit. Values are not transferable to systems not listed on the assay sheet. The product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product.

CELL-DYN Enzymatic Cleaner Concentrate solution must be thoroughly removed from the analyzer Wash Block and/ or Probe before running control products.

Control products are not to be used as calibrators. The white blood cell components simulate white blood cells in size, not morphology.

EXPECTED VALUES AND THEIR DERIVATION

Refer to the table of values on the enclosed assay sheet for expected results. The assay values are obtained from replicate testing on CELL-DYN Systems. The mean range is an estimate of observed interlaboratory variation due to reagent differences, maintenance, calibration, and operating technique. The Systems are calibrated as follows:

- using whole blood according to CLSI and ICSH-recommended manual, reference methods^{5,6,7,8,9,10}
- using a commercially available CELL-DYN calibrator
- using the appropriate, reliably calibrated CELL-DYN hematology analyzer.

PERFORMANCE CHARACTERISTICS

The assay values, with their associated mean ranges, reflect the expected biological variability of the control materials and the estimated interlaboratory variation.

Each laboratory should establish a mean and acceptable range for each lot of control material. The laboratory mean should fall within the listed Mean Range. An individual laboratory's range may include results above and below the listed Mean Range. Laboratories may consider results acceptable when at least 95% of test results are within the laboratory's expected range. For additional information, refer to the appropriate CELL-DYN Operator's Manual.

BIBLIOGRAPHY

- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. Laboratory Biosafety Manual. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- The Expert Panel in Cytometry on the ICSH, J.M. England, et.al.: Recommended Methods for the Visual Determination of White Cell and Platelet Counts. World Health Organization, WHO/LAB/88.3, 1988.
- The Expert Panel in Cytometry on the ICSH, J.M. England, et.al.: Reference Method for the Enumeration of Erythrocytes and Leucocytes. Clin. Lab. Haemat. 16:131-138, 1994.
- Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard - Third Edition; CLSI Document H15-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2000.
- Procedure for Determining Packed Cell Volume by the Microhematocrit Method; Approved Standard - Third Edition; CLSI Document H07-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2000.
- The Expert Panel in Cytometry on the ICSH, J.M. England, et.al.: The Assignment of Values to Fresh Blood Used for Calibrating Automated Blood Cell Counters. Clin. Lab. Haemat. 10:203-212, 1987.
- Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard - Second Edition; CLSI Document H26-A2, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2010.

ORDERING INFORMATION

[REF] 08H58-01	CELL-DYN 29 Plus Control (with Retic)	12 x 3.0-mL Vials
[REF] 08H58-02	CELL-DYN 29 Plus Control (with Retic)	6 x 3.0-mL Vials

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BIBLIOGRAFIA

- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- The Expert Panel in Cytometry on the ICSH, J.M. England, et.al.: Recommended Methods for the Visual Determination of White Cell and Platelet Counts. World Health Organization, WHO/LAB/88.3, 1988.
- The Expert Panel in Cytometry on the ICSH, J.M. England, et.al.: Reference Method for the Enumeration of Erythrocytes and Leucocytes. Clin. Lab. Haemat. 16:131-138, 1994.
- Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard - Third Edition*; CLSI Document H15-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2000.
- Procedure for Determining Packed Cell Volume by the Microhematocrit Method; Approved Standard - Third Edition*; CLSI Document H07-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2000.
- The Expert Panel in Cytometry on the ICSH, J.M. England, et.al.: The Assignment of Values to Fresh Blood Used for Calibrating Automated Blood Cell Counters. Clin. Lab. Haemat. 10:203-212, 1987.
- Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard - Second Edition*; CLSI Document H26-A2, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2010

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REF 08H58-01	CELL-DYN 29 Plus Control (with Retic)	12 provette da 3,0 mL
REF 08H58-02	CELL-DYN 29 Plus Control (with Retic)	6 provette da 3,0 mL

Key to symbols used/
Erläuterung der verwendeten Symbole/
Légende des symboles utilisés/
Símbolos utilizados/
Legenda dei simboli utilizzati

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Abbott Laboratories

Abbott

© 2005, 2018 Abbott Laboratories. Printed in the USA.
9231565C 350490-11 August 2018

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- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.
- The Expert Panel in Cytometry on the ICSH, J.M. England, et.al.: Recommended Methods for the Visual Determination of White Cell and Platelet Counts. World Health Organization, WHO/LAB/88.3, 1988.
- The Expert Panel in Cytometry on the ICSH, J.M. England, et.al.: Reference Method for the Enumeration of Erythrocytes and Leucocytes. Clin. Lab. Haemat. 16:131-138, 1994.
- Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard - Third Edition*; CLSI Document H15-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2000.
- Procedure for Determining Packed Cell Volume by the Microhematocrit Method; Approved Standard - Third Edition*; CLSI Document H07-A3, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2000.
- The Expert Panel in Cytometry on the ICSH, J.M. England, et.al.: The Assignment of Values to Fresh Blood Used for Calibrating Automated Blood Cell Counters. Clin. Lab. Haemat. 10:203-212, 1987.
- Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard - Second Edition*; CLSI Document H26-A2, CLSI, 940 West Valley Rd., Suite 1400, Wayne, PA 19087, USA, 2010.

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