

CONTROL L N H

 Exp. 2024-04-05

SYSTEM	CONTROL L		CONTROL N		CONTROL H	
	LOT	L40229	LOT	N40229	LOT	H40229
PARAMETER	TARGET	± RANGE **	TARGET	± RANGE **	TARGET	± RANGE **
WBC 10 ⁹ /L	2.94	0.56	6.63	1.26	15.8	3.1
NEU 10 ⁹ /L	1.30	0.27	3.19	0.64	8.11	1.60
%N	44.3	6.0	48.1	5.0	51.4	6.0
LYM 10 ⁹ /L	0.89	0.20	1.70	0.50	3.45	1.00
%L	30.1	8.0	25.7	6.0	21.9	5.0
MONO 10 ⁹ /L	0.34	0.20	0.73	0.40	1.69	0.60
%M	11.4	6.0	11.1	5.0	10.7	3.0
EOS 10 ⁹ /L	0.10	0.10	0.19	0.17	0.44	0.19
%E	3.00	3.00	2.80	2.00	2.80	1.00
BASO 10 ⁹ /L	0.10	0.10	0.25	0.25	0.50	0.50
%B	1.50	1.50	1.50	1.50	1.50	1.50
RBC 10 ¹² /L	2.73	0.25	4.05	0.39	5.03	0.83
HGB g/dL	7.47	0.30	11.5	0.5	16.3	0.8
HCT %	22.7	3.4	34.4	5.4	46.8	11.8
MCV fL	83.2	7.0	84.9	7.5	93.0	9.2
MCH pg	27.4	2.6	28.4	2.6	32.3	4.9
MCHC g/dL	32.9	4.9	33.4	5.0	34.8	7.8
RDW %	15.9	2.5	16.2	2.5	15.8	2.5
NRBC 10 ⁹ /L*	0.01	0.01	0.01	0.01	2.20	0.60
NR/W*	0.20	0.20	0.08	0.08	14.0	2.5
PLT 10 ⁹ /L	66.6	20.0	210	50	459	62
MPV fL	6.96	2.00	6.91	2.00	6.91	2.00
RETIC 10 ⁹ /L	224	50	162	50	143	50
%R	8.20	1.50	4.01	1.00	2.84	0.80
IRF	0.60	0.18	0.49	0.14	0.35	0.10
IG 10 ⁹ /L	0.31	0.11	0.76	0.22	1.95	0.43
%IG	10.5	3.4	11.5	3.0	12.4	2.2
%rP	9.92	2.00	9.02	1.70	9.27	1.60
MCHr pg	27.4	2.7	29.9	3.3	35.8	4.2

SYSTEM	CONTROL L		CONTROL N		CONTROL H	
	LOT	L40229	LOT	N40229	LOT	H40229
PARAMETER	TARGET	± RANGE **	TARGET	± RANGE **	TARGET	± RANGE **
WBC 10 ⁹ /L	3.04	0.96	6.78	1.04	16.1	2.0
NEU 10 ⁹ /L	1.34	0.60	3.26	0.85	8.29	1.80
%N	44.0	6.2	48.1	5.3	51.6	5.0
LYM 10 ⁹ /L	0.93	0.47	1.73	0.53	3.47	0.80
%L	30.7	4.8	25.6	3.6	21.6	2.2
MONO 10 ⁹ /L	0.35	0.25	0.77	0.40	1.77	0.70
%M	11.4	4.0	11.4	4.1	11.0	2.8
EOS 10 ⁹ /L	0.10	0.10	0.19	0.15	0.45	0.18
%E	3.00	3.00	2.81	1.16	2.80	0.71
BASO 10 ⁹ /L	0.10	0.10	0.25	0.25	0.50	0.50
%B	1.50	1.50	1.50	1.50	1.50	1.50
RBC 10 ¹² /L	2.80	0.34	4.13	0.43	5.11	0.70
HGB g/dL	7.60	0.54	11.7	0.6	16.4	0.7
HCT %	24.5	3.4	36.9	5.1	49.9	8.7
MCV fL	87.6	5.6	89.4	6.2	97.8	6.6
MCH pg	27.1	4.1	28.3	3.4	32.2	4.8
MCHC g/dL	31.0	7.2	31.6	6.4	32.9	7.9
RDW %	12.2	2.5	12.4	2.5	11.5	2.5
NRBC 10 ⁹ /L*	0.01	0.01	0.01	0.01	2.25	0.50
NR/W*	0.20	0.20	0.08	0.08	14.0	1.8
PLT 10 ⁹ /L	68.2	30.0	210	24	458	66
MPV fL	9.12	2.00	9.16	2.00	9.14	2.00
RETIC 10 ⁹ /L	227	50	162	45	133	58
%R	8.12	1.50	3.92	1.00	2.61	0.80
IRF	0.33	0.18	0.28	0.14	0.21	0.10
IG 10 ⁹ /L	0.31	0.23	0.76	0.35	1.94	0.91
%IG	10.3	4.1	11.2	3.3	12.1	4.1
%rP	10.5	2.2	8.96	0.90	9.04	0.60
MCHr pg	23.8	2.7	26.8	3.3	32.4	4.2

